

### **Embryo Fixation Protocol:**

- 1) Collect embryos on grape juice agar plate and age to appropriate stage.
- 2) Using a paint brush, move embryos from the collection plate to a basket in embryo wash solution (7% NaCl, 0.7% Triton X-100).
- 3) Dechorionate by incubated basket of embryos in 50% bleach for 2 min.
- 4) Wash embryos in basket again.
- 5) With a squirt bottle, wash embryos down to one side of the basket, and blot dry with a kimwipe from underneath.
- 6) Transfer embryos to an eppendorf (or other tube) by squirting heptane on the embryos from the bottom of the basket.
- 7) Add an equal volume of fix solution to the heptane. You want 1:1 heptane:fix and your embryos floating at the interface.
- 8) Fix by rocking on a Nutator for the appropriate time (see below).
- 9) Remove the fix by pipetting or aspiration, leaving the heptane behind. Make sure you remove all or most of the fix, otherwise they won't "pop" at the next step.
- 10) Add an equal volume of MeOH (1:1 MeOH:Heptane) and immediately shake vigorously. The devittelinized embryos will sink to the bottom, and the ones that don't pop will remain at the interface.
- 11) Aspirate off everything except the embryos at the bottom.
- 12) Wash 2 X with MeOH and 2 X with EtOH, and store at -20°C in EtOH.

FIX: straight 37% formaldehyde for 5 min

OR

7% formaldehyde in PBS for 15-25 min depending on application.