

## EMBRYO IN SITU HYBRIDIZATION PROTOCOL USING DIG-RNA PROBES

1. Collect, dechorionate, fix, and devittelinize embryos  
Fix either 37% (straight) formaldehyde:heptane for 5 min  
or 7% formaldehyde in PBT:heptane for 25 min. Fixation is important to  
detect mRNA, and you may have to optimize yourself.
2. You can store embryos in ethanol at -20°C indefinitely.
3. When ready to prehybe, transfer embryos to MeOH
4. Rinse with 1:1 MeOH:5% formaldehyde/PBT
5. Rinse with 5% formaldehyde/PBT
6. Post-fix 15 min with 5% formaldehyde/PBT
7. Wash 5 X 5 min with PBT  
important to wash all 5 times and for 5 min; need to get rid of fixative
8. Wash 10 min in 1:1 PBT:hybe buffer
9. Wash 10 min in 1:1 hybe buffer
10. Prehybe at least 1 hour in hybe buffer at 70°C  
Use about 100ul for 30ul volume of embryos  
No need to rock, just put eppendorf in a temp block with water.
11. Add Dig-RNA probe in hybe to embryos  
Again, use about 100ul; just make sure embryos are well covered
12. Hybe overnight at 70°C  
Again, just in the temp block or water bath
13. In morning remove and discard probe and wash:  
1 X 45 min hybe  
1 X 45 min 66:33 hybe:PBT  
1 X 45 min 33:66 hybe:PBT  
4 X 20 min PBT  
All these washes at 70°C, and warm up solutions first; no need to agitate;  
just use the temp block or water bath
14. Incubate embryos in 1:2000 anti-Dig alkaline phosphatase in PBT  
anti-Dig AP from Boehringer Mannheim; no need to preadsorb antibody or  
to block the embryos; rock on Nutator
15. Wash 3 X 20 min with PBT on Nutator  
I've also used 10 min washes

16. Rinse 3 X 2 min with reaction solution:
  - 100 mM Tris pH 9.5
  - 100 mM NaCl
  - 50 mM MgCl<sub>2</sub>
  - 0.1% Tween 20
17. Color reaction:
  - Add 4.5ul NBT and 3.5ul X phosphate per ml of reaction solution
  - Incubate embryos in this in the dark for 5 min to overnight (often O/N is needed)
18. Stop reaction by washing several times with PBT

Hybe solution:

- 50% deionized formamide
- 5X SSC
- 100 ug/ml sonicated, denatured salmon sperm DNA
- 100 ug/ml E. coli tRNA
- 50 ug/ml heparin
- 0.1% Tween 20
- pH to 4.5 with Citric Acid

NBT= 4-nitro blue tetrazolium chloride  
X-phosphate= 5-bromo-4-choro-3-indolyl-phosphate  
Buy as solutions from Boehringer and store -20°C

Transcribe RNA probe using B.M. Dig labeling kit  
No need to fragment the probe, although some people say it helps  
Correct concentration of probe found by doing dilution series on embryos  
Usually 1:2000 works for me

For example:

- 2ul=1ug cut plasmid DNA
- 2ul 10 X reaction buffer (Boehringer)
- 2ul 10 X dig-UTP + NTP mix (Boehringer)
- 1ul RNAsin (Boehringer)
- 1ul T7 or T3 RNA pol (Boehringer)
- up to 20ul with water

2 hours at 37°C  
run 1 ul of reaction on a gel

Steps 10-13 (70°C ones) no need to agitate embryos. Just let them sit in a temp block or water bath in an eppendorf tube. I usually rock the embryos on a nutator through all the other steps. Hybe temp is important: make sure the embryos are really at 70°C. If this doesn't work for your probe, 65°C might be better.