

CELL SORTING AND HOECHST STAINING ON ISOLATED FOLLICLE CELLS. Pelin Cayirlioglu/Duronio lab
Based on Bryant et al PNAS 96:5559-5564 (1999).

- 1) Dissect out 80-150 females (fat) in Schneider's medium (you don't need to fine dissect them, just get the two ovaries out of the females).
- 2) Transfer them into an eppendorf tube in Schneider's medium.
- 3) Remove the Schneider's and wash the ovaries with 1XPBS once.
- 4) Incubate the ovaries in 0.5% trypsin (1:250 sigma) for 15 minutes, with intermittent vortexing.
- 5) Filter the supernatant (that has the dissociated follicle cells) through a 80 μ mesh into 1ml of Schneider's solution.
- 6) Spin the supernatant at 4200rpm for 7 minutes to pellet the cells.
- 7) Resuspend the pellet in 200-300 μ l Schneider's and keep on ice.
- 8) Repeat steps 4-7 2-3 more times, until you stop getting any cells in the supernatant (you will understand this because the supernatant will no longer be blurry if there are no cells in there).
- 9) After you got all the cells, if you want to do some FACS on them, stain them with 5 μ g/ml Hoechst 33342 for 45 minutes on the nutator at room temperature. Store on ice until sorting.
- 10) If you want to get RNA/DNA from them, spin down all the cells, and you can start the extraction step immediately or you can save the cell pellet at -80.