

## X-GAL Staining of Embryos

To obtain homozygous mutant embryos for PCR analysis. (Make sure your stock is balanced with a lacZ balancer chromosome; e.g. wg-lacZ, Ubx-lacZ, ftz-lacZ.)

- 1) Collect germ band retracted embryos (i.e. older than 8 hours at 25°C). Overnight collections also work, but pick out germ band retracted (i.e. older) embryos.
- 2) Fix the embryos in one of two ways:
  - a. with MeOH/Heptane (i.e. just pop the embryos with MeOH)
  - b. with 1% glutaraldehyde/PBS 1:1 with heptane for 15 minutes
- 3) Devitellinize (pop) the embryos with MeOH.
- 4) Wash 3 X with 1 ml of MeOH (just add the MeOH and let embryos sink).
- 5) Wash 3 X 1 min with 1 ml PBT.
- 6) Rock embryos on nutator in 1 ml of PBT for 6 hours (or until the end of the day).
- 7) Wash once with 1 ml staining solution, and then rock O/N in 1 ml staining solution in the dark (i.e. wrap the tubes in foil).

All steps at room temperature

X-GAL staining solution for 2 ml:    1780 ul PBT  
    150 ul 2% X-GAL in DMF  
    70 ul Ferri-ferro cyanide  
    2 ul 1M MgCl<sub>2</sub>

Ferri-ferro cyanide:            0.1M K ferricyanide  
    0.1M K ferrocyanide

Keep at 4°C in the dark (we have some in the refrigerator).

- 8) After staining wash 3 X with PBT
- 9) Hand pick the white embryos (5-10) with a small pipet man, and put in 0.5 ml tube.
- 10) Remove all the PBT, replace with 20 ul squish solution, and mash with rounded pipet tip.

Squish solution:            10 mM Tris 8.0  
    1 mM EDTA  
    25 mM NaCl  
    200 ug/ml Proteinase K
- 11) Heat the mashed embryos at 37°C for 15 min
- 12) Inactive the proteinase K by heating at 85°C for 15 min
- 13) Spin down junk. Use 1-5 ul for PCR